COMPOSITION

Hyaluronic Acid gel, 45 mg/mL - 22.5 mg (H-HA) + 22.5 mg (L-HA)/mL
Phosphate buffered saline, pH 7.0 q.s.
One syringe contains 2.0mL of Neurafill® Revive

DESCRIPTION

Neurafill® Revive is a sterile, transparent, pyrogen-free, viscoelastic gel of non-animal origin hyaluronic acid. This product is presented in a pre-filled, disposable glass syringe with a volume of 2.0 mL for single use only.

Each box of Neurafill® Revive contains one 2.0 mL pre-filled syringe, one sterile (22G x 40 mm) cannula and a set of labels including batch number, one of which should be attached to the patient's file and the other should be given to the patient in order to ensure traceability.

STERILIZATION

The content of the Neurafill® Revive syringes is sterilized by moist heat. The provided needles and cannulas are sterilized by radiation.

INDICATIONS

Neurafill® Revive is intended for facial treatment, restoring superficial fat compartment

CONTRAINDICATIONS

- Neurafill® Revive must not be used in concurrence with treatments such as laser resurfacing and medium deep skin-peeling.
- Do not inject intravenously, into muscles, tendons or for mammary augmentation.
- Neurafill® Revive must not be used in:
 - Patients who tend to develop hypertrophic scarring and have keloid-prone skin
 - Patients who are known to be hypersensitive to hyaluronic acid
 - Patients who have a history of autoimmune disease
 - Patients who have coagulation disorders
 - Patients who have streptococcus skin infections
 - Women who are pregnant or breast-feeding
 - Patients under 18 years old
- Neurafill® Revive must not be injected into areas with inflammatory and/or infectious problems.

PRECAUTIONS FOR USE

- Do not inject Neurafill® Revive into an area which has already been treated with filling product within three (3) months.
- Neurafill® Revive must be injected into healthy, noninflamed and disinfected skin.
- Patients who are showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a dual test before any injection is administered. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.
- Patients who are undergoing anti-coagulant treatment must be warned of the increased risk of hematomas and bleeding during the injection.
- It is recommended to avoid taking aspirin or high doses of vitamin E few days before the injection.
- Advise the patients to avoid UV exposure for 3-5 days after the injection. SPF50 sunscreen cream can be used starting second day.
- If the cannula is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the cannula.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. Therefore, Neurafill® Revive should never be placed in contact with such substances or with medical-surgical instrumentation which has been treated with this type of substance.

UNDESIRABLE EFFECTS


Before treatment, medical practitioners must inform the patient that potential side effects associated with use of Neurafill® Revive can occur immediately or may be delayed.

These include:

- Inflammatory reactions (redness, edema, erythema, etc.) that may be associated to itching, pain due to pressure can occur following injection. These reactions may last few hours, days and in some cases for a week. These symptoms can be relieved by applying ice to the treated area. They generally disappear in a short period of time but may persist for up to 10 days.
- In extremely rare cases, as described in the literature discoloration at the injection site, necrosis of the glabellar area (between the eyebrows above the nose), abscess, granuloma and hypersensitivity may occur after injections, so injection should be performed by professional doctor.
- Local bleeding or hematomas in the injected area.
- Patients must report all secondary reactions which persist for more than one week to medical practitioner as soon as possible. The medical practitioner should use an appropriate treatment.
- Any other undesirable side effects associated with injection of Neurafill® Revive must be reported to the distributor and/or the manufacturer.

METHOD OF USE – POSOLOGY

- Neurafill® Revive can be applied through the most common injection techniques using cannula in the superficial fat layer of the skin.
- Before beginning treatment, patient must be informed of product's indications, contra-indications, incompatibilities and potential undesirable effects.
- The area to be treated must be disinfected properly before the injection.
- Massage gently around the injection site.
- An initial cycle of two treatment sessions at 30-day intervals is recommended, followed if necessary by maintenance treatments every 2 months using by Neurafill® Aqua or Neurafill® Revive. However, it is suggested to evaluate the specific Neurafill® Revive protocol according to the patient's degree of aging.
- It is very important that the cannula is properly assembled to the syringe. Improper assembly may cause separation of the cannula and syringe during injection.
- Grasp Luer lock tip cap, break it (Fig. 1).
- Screw the cannula clockwise until it is fully locked (Fig. 2).
- Remove the protective cap by pulling it straight off the syringe (Fig. 3).
- Inject Neurafill® Revive at room temperature using strict aseptic methods.
- If the cannula becomes obstructed and the injection pressure become too high, stop the injection and change the cannula.
- After the injection, the medical practitioner should massage the treated area gently in order to distribute the product uniformly.

WARNING

- Do not use if package is damaged
- Check the expiry date on the product label
- Do not mix with other products
- Discard the syringe and remaining product after use
- Do not re-sterilize
- Do not re-use

STORAGE CONDITIONS

- The expiry date is indicated on package
- Store between 2°C and 25°C
- Protect from freezing, heat and sunlight
- Fragile

Manufacturer and Authorization Holder:
ARYA TAJHIZ IBN SINA
Number 21, 20 St, Kaveh Industrial City, Saveh, Iran

SYMBOLS

	Do not reuse
	Caution
	Sterilized using steam
	Storage temperature (2 ~ 25°C)
	Batch code
	Use-by date
	Do not use if package is damaged
	Consult instructions for use
	Do not re-sterilize
	Keep away from sunlight
	Date of manufacture
	Manufacturer

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