Neurafill® Lidocaine
Instructions for Use

EN

COMPOSITION
Hyaluronic Acid gel 20 mg/mL
Lidocaine hydrochloride 3 mg/mL
Phosphate buffered saline, pH 7.0 q.s.
One syringe contains 1.0 mL of Neurafill® Lidocaine

DESCRIPTION
Neurafill® Lidocaine is a sterile, transparent, pyrogen-free, viscoelastic gel of cross-linked non-animal origin hyaluronic acid. This product is presented in a pre-filled, disposable glass syringe with a volume of 1.0 mL for single use only.
Each box of Neurafill® Lidocaine contains one 1.0 mL pre-filled syringe, two sterile 30G 1/2" disposable needles and a set of labels including batch number, one of which should be attached to the patient's file and the other should be given to the patient in order to ensure traceability.

STERILIZATION
The content of the Neurafill® Lidocaine syringes is sterilized by moist heat. The provided needles are sterilized by radiation.

INDICATIONS
Neurafill® Lidocaine is designed to be injected into the dermis for the treatment of wrinkles and folds on the lower face, specifically the nasolabial folds.
The presence of lidocaine is meant to reduce the patient's pain during treatment.

CONTRA-INDICATIONS
• Do not inject Neurafill® Lidocaine into the eye circles and eyelids.
• Do not inject into the blood vessels.
• Do not overcorrect.
• Neurafill® Lidocaine must not be used in:
- Patients who tend to develop hypertrophic scarring
- Patients with known hypersensitivity to lidocaine or to amide type local anesthetics
- Patients who have a history of autoimmune diseases
- Patients who have coagulation disorders
- Patients who have streptococcus linked dermal disease
- Women who are pregnant or breast-feeding
- Children
• Neurafill® Lidocaine must not be injected into areas with inflammatory and/or infectious problems.
• Neurafill® Lidocaine must not be used in association with a treatment such as laser therapy, chemical peeling or dermabrasion.

PRECAUTIONS FOR USE
• Do not use Neurafill® Lidocaine for injections other than intradermal injections.
• Do not inject Neurafill® Lidocaine into an area which has already been treated with another permanent filling product.
• Neurafill® Lidocaine must be injected into healthy, noninflamed and disinfected skin.
• Patients who are showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a dual test before any injection is administered. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.
• Patients who are undergoing anti-coagulant treatment must be warned of the increased risk of haematomas and bleeding during the injection.
• It is recommended to avoid taking aspirin or high doses of vitamin C the week before the injection.
• It is recommended not to apply make up for 24 hours after injection and to avoid prolonged exposure to sunlight, UV light, cold temperatures or using saunas or steam baths for the two weeks after injection.
• If the needle is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the needle.
• Athletes should be made aware that this product contains an active principle that may produce a positive result in antidoping tests.
• Medical practitioners must take into account the fact that this product contains lidocaine.
• Considerations should be given to the total dose of lidocaine administered if dental block or topical administration of lidocaine is used concurrently. High dose of lidocaine (more than 400 mg) can cause acute toxic reactions manifesting as symptom affecting the central nervous system and cardiac conduction.

INCOMPATIBILITIES
Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. Therefore, Neurafill® Lidocaine should never be placed in contact with such substances or with medical-surgical instrumentation which has been treated with this type of substance.

UNDESIRABLE EFFECTS
Before treatment, medical practitioners must inform the patient that potential side effects associated with use of Neurafill® Lidocaine can occur immediately or may be delayed.
These include:
• Inflammatory reactions (redness, edema, erythema, etc.) that may be associated to itching, pain due to pressure can occur following injection. These reactions may last few hours, days and in some cases for a week.
• Local bleeding or haematomas in the injected area.
• Induration, staining and poor effect at the injection site.
• Cases of necroses in the glabellar region, granulomas, hypersensitivity and abscess formation have been reported in the other products after the injection of hyaluronic acid.
• Patients must report all secondary reactions which persist for more than one week to medical practitioner as soon as possible. The medical practitioner should use an appropriate treatment.

• Any other undesirable side effects associated with injection of Neurafill® Lidocaine must be reported to the distributor and/ or the manufacturer.
The following undesirable effects associated with lidocaine injections have been reported:
• Shock: Medical practitioners must do close observation for occurring the shock. If Patients showing decreased blood pressure, pallor of face, abnormality of pulse, respiratory depression, etc., the medical practitioner should discontinue the injection as soon as possible and use an appropriate treatment.
• Malignant hyperthermia: Severe malignant hyperthermia accompanied by cryptogenic tachycardia, arrhythmia, blood pressure fluctuation, rapid increase of body temperature, muscle rigidity, cyanosis, hyperventilation, diaphoresis, acidosis and/or hyperkalemia, myoglobinuria (erythruia) are may occur rarely. Should be stopped immediately if these symptoms are appeared when Neurafill® Lidocaine injection and then use an appropriate treatment. (Intravenous injection of dantrolene sodium, general hypothermy, hyperventilation by pure oxygen, correction of acid-base equilibrium, etc.). In this case, urine output should be monitored because these symptoms may lead the renal failure.
• Central nervous system:
- Medical practitioners must do close observation for occurring the tremor, convulsion, etc. Should be stopped the injection immediately if these symptoms are appeared and then use an appropriate treatment (Diazepam, ultrashortacting barbiturate (thiopental sodium), etc.).
- Medical practitioners must do close observation for occurring the somnolence, anxiety, irritability, inattention, dizziness, nausea, vomiting, etc. and then use an appropriate treatment if necessary.
• Anaphylaxis: Skin symptoms such as urticaria or swelling are may occur.

METHOD OF USE – POSOLOGY
• Neurafill® Lidocaine is designed to be injected into the mid dermis to deep dermis for the treatment of wrinkles and folds on the lower face, specifically the nasolabial folds by medical practitioners who have undertaken specific training in the injection technique for filling.
• Before beginning treatment, patient must be informed of product's indications, contra-indications, incompatibilities and potential undesirable effects.
• The correct injection technique of Neurafill® Lidocaine is important for the effectiveness.
- The area to be treated must be disinfected properly before the injection.
- It is recommended to perform injections using the needles provided in the box.
- To avoid needle breakage do not attempt to bend the needle.
- Neurafill® Lidocaine should be injected into the dermis
- A too large volume or a too superficial injection may give bumps on the treatment site and if any bumps or irregularities are noted after injections, the medical practitioner should be gently massaged the treated area immediately.
• It is very important that the needle is properly assembled to the syringe. Improper assembly may cause separation of the needle and syringe during injection.
- Grasp Luer lock tip cap and unscrew it (Fig. 1).
- Screw the needle clockwise until it is fully locked (Fig. 2).
- Remove the protective cap by pulling it straight off the syringe (Fig. 3).
- Start the injection slowly into the mid dermis.
• If the needle becomes obstructed and the injection pressure become too high, stop the injection and change the needle.
• After the injection, the medical practitioner should massage the treated area gently in order to distribute the product uniformly.

Fig 1. Grasp Luer lock tip cap and unscrew it.


Fig 2. Screw the needle clockwise until it is fully locked (Grasp the lure lock adapter of the syringe barrel firmly. Press and turn the needle shield clockwise until the needle locks onto the syringe).



Fig 3. Remove the protective cap by pulling it straight off the syringe

WARNING
• Do not use if package is damaged
• Check the expiry date on the product label
• Do not mix with other products
• Discard the syringe and remaining product after use
• Do not re-sterilize
• Do not re-use

STORAGE CONDITIONS
• The expiry date is indicated on package
• Store between 2°C and 25°C
• Protect from freezing, heat and sunlight
• Fragile

Manufacturer and Authorization Holder:
ARYA TAJHIZ IBN SINA
Number 21, 20 St, Kaveh Industrial City, Saveh, Iran

Table with 2 columns: Symbol and Description. Symbols include: No reuse, Caution, Sterile, Storage temperature, Batch code, Use-by date, Do not use if package is damaged, Consult instructions for use, Do not resterilize, Keep away from sunlight, Date of manufacture, Manufacturer.

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